

Quick, Effective, Non-Invasive Skin Closure

- **Fast & Strong:** Close high tension wound in seconds, yet easy removal with no skin damage
- **Simple & Safe:** Mastered in minutes, with no needles, puncture pain, burning or lidocaine
- **Multi-Functional:** From the frontline to the hospital; for temporary closure or final closure
- **Versatile:** Usable on straight, jagged or curving wounds; range of sizes and adhesives
- **Usable on Wet Wounds:** Can be applied to wet wounds & eversion helps with hemostasis
- **Tactile:** Safely apply in low visibility with tactile feedback
- **Force Multiplier:** Helps keep warfighter on the battlefield
- **72-Hour Sustainability:** Helps sustain the wounded for 72+ hours, lasts 5+ days
- **All Roles of Care:** Works for wound care in ROC 1, ROC 2, and ROC 3
- **Lightweight & Portable:** Compact and nearly flat; light - only 0.15 oz

Simple Process: “Pick, Place & Pull”



Unique Entity ID: VC6SGH5HDBL5

NAICS Code: 339113, Surgical Appliance & Supplies Manufacturing

CAGE Code: 7K3F7

NSNs: (See chart on reverse side)

DAPA & ECAT: Beacon Point Associates, LLC (SDVOSB)
 DAPA: SPO200-18-H-0042
 ECAT: SPE2DE-21-D-7004

Contact: DermaClip Team, Military & Government
 MIL@dermaclip.com
 Beacon Point Associates, LLC (SDVOSB)
 Quotes@TeamBPA.com / +1 (239) 673-6965

Headquarters: DermaClip US, LLC
 730 N. Post Oak Road, Ste. 202
 Houston, TX 77024
 +1 (713) 682-3185 / +1 (844) 572-7200
 www.DermaClip.com



DermaClip® gives providers options related to clip size and target wound length to allow customized closure for each wound or incision. **Please note the Extreme product utilizes an extra strength adhesive not suitable for fragile skin.**

DermaCLIP® Sizes & Packaging Setups

Clip Size	Target Closure Length 1 INCH / PKG		Target Closure Length 3 INCHES / PKG	
	REGULAR CLIP	LARGE CLIP	REGULAR CLIP	LARGE CLIP
	62001	62003	62002	62004
Box Catalog No. (12 Pkgs) Case Catalog No. (144 Pkgs)	63001	63003	63002	63004
Description	Package Contents		Box 12 Pkg	Case 12 Boxes
Regular Clip (11 mm) 2 Clips / Package Standard Adhesive			MPN 62001 NSN PENDING	63001 PENDING
Regular Clip (11 mm) 5 Clips / Package Standard Adhesive			MPN 62002 NSN 6515-01-662-6938	63002 6515-01-662-6997
Large Clip (21 mm) 1 Clip / Package Standard Adhesive			MPN 62003 NSN PENDING	63003 PENDING
Large Clip (21 mm) 3 Clips / Package Standard Adhesive			MPN 62004 NSN 6515-01-661-8487	63004 6515-01-662-6961
Large Clip (21 mm) 3 Clips / Package EXTREME™ Adhesive			MPN 62005 NSN 6515-01-673-7948	63005 6515-01-673-7942

DAPA Beacon Point Associates, LLC, an SDVOSB (Contract # SPO200-18-H-0042)

ECAT Beacon Point Associates, LLC, an SDVOSB (Contract # SPE2DE-21-D-7004)

NSN (chart) **Unique Entity ID** VC6SGH5HDBL5 **CAGE** 7K3F7 **NAICS** 339113, Surgical Appliance & Supplies Mfg.

All products are Rx Only devices and are EtO sterilized. Standard Adhesive products have a 3-year shelf life & EXTREME™ Adhesive have a 4-year shelf life from date of manufacture.