


ONE-HANDED CLOSURE FOR THE WARFIGHTER

QUICK, EFFECTIVE, NEEDLE-FREE WOUND CLOSURE

- **Easy to Use:** Just place, press together & pop up to close.
- **Fast & Strong:** Close minor wound in seconds.
- **Easy Removal:** No skin damage, good for transport use.
- **Tactile:** Easy to apply in low visibility with tactile feedback.
- **Force Multiplier:** Warfighter can apply with no assistance.
- **Long Duration & Durable:** Stays on for 5+ days; can get wet.
- **ROC 1 Care:** Works for self-aid, buddy-aid, TCC, MEDEVAC, BAS.
- **Lightweight & Portable:** Nearly flat; light - only 0.15 oz.



THE WARFIGHTER'S SOLUTION TO MINOR WOUNDS



**DERMAQUIK
OTC CLIP**

**PLACE IT
ON THE SKIN**

**PRESS IT
TOGETHER**

**POP UP
TO CLOSE**

Description		Package Contents	Box 12 Pkg	
	OTC, 1 Clip/Pkg, Standard Adhesive		MPN 62501 NSN PENDING	
	OTC, 1 Clip/Pkg, EXTREME Adhesive		MPN 62502 NSN PENDING	

DAPA Beacon Point Associates, LLC, an SDVOSB (Contract # SPO200-18-H-0042)

ECAT Beacon Point Associates, LLC, an SDVOSB (Contract # SPE2DE-21-D-7004)

NSN (chart) **Unique Entity ID** VC6SGH5HDBL5 **CAGE** 7K3F7 **NAICS** 339113, Surgical Appliance & Supplies Mfg.

MFR. Contact: DermaQuik Team, Military & Gov't. MIL@dermaquik.com. +1 (713) 682-3185 / +1 (844) 572-7200

MFR. Headquarters: DermaClip US, LLC. 730 N. Post Oak Rd., Ste. 202, Houston, TX, 77024, USA

Contract Holder: Beacon Point Associates, LLC (SDVOSB). Quotes@TeamBPA.com / +1 (239) 673-6965

All products are EtO sterilized. Home & Field products have a 3-year shelf life. Sport & Extreme products have a 4-year shelf life from date of manufacture. For more info, visit www.DermaQuik.com.